



Arthroscopic instrument
关节镜器械



QUALITY
OUR WAY TO SUCCESS

Arthroscopic instrument

篮钳
Basket Forceps



	L × B	鸭嘴, 直	straight
JP2003RB	125 × 3.2,	鸭嘴, 直	straight
JP2007RB	125 × 3.2,	鸭嘴, 上翘	upwards angled
JP2008RB	125 × 3.2,	鸭嘴, 左弯, 上翘	left curved, upwards
JP2009RB	125 × 3.2,	鸭嘴, 右弯, 上翘	right curved, upwards
JP2100RB	82 × 2.3,	鸭嘴, 直	straight
JP2101RB	82 × 2.3,	鸭嘴, 上翘	upwards angled
JP2120RB	215 × 3.2,	鸭嘴, 直	straight
JP2121RB	215 × 3.2,	鸭嘴, 上翘	upwards angled
JP2124RB	215 × 3.2,	鸭嘴, 上弯	upwards curved



关节镜器械
Arthroscopic instrument



JP2020RB	125 × 3.2,	鸭嘴, 直	straight
JP2021RB	125 × 3.2,	鸭嘴, 左弯	left curved
JP2022RB	125 × 3.2,	鸭嘴, 右弯	right curved
JP2023RB	125 × 3.2,	鸭嘴, 上弯	upwards curved



JP2050RB	125 × 1.8,	窄形, 上翘	upwards angled
JP2058RB	125 × 1.8,	窄形, 直	straight
JP2059RB	125 × 1.8,	窄形, 上翘	upwards angled
JP2110RB	82 × 1.8,	窄形, 直	straight
JP2127RB	215 × 1.8,	窄形, 直	straight
JP2128RB	215 × 1.8,	窄形, 上翘	upwards angled

JP2038RB	125 × 4.4,	卵圆形, 直	straight
JP2039RB	125 × 4.4,	卵圆形, 左弯	left curved
JP2040RB	125 × 4.4,	卵圆形, 右弯	right curved
JP2042RB	125 × 4.4,	卵圆形, 上翘	upwards angled
JP2105RB	82 × 2.3,	卵圆形, 直	straight
JP2106RB	82 × 2.3,	卵圆形, 左弯	left curved
JP2107RB	82 × 2.3,	卵圆形, 右弯	right curved



JP2031RB	125 × 1.8,	窄形, 直	straight
JP2071RB	125 × 2.2,	窄形, 上翘	upwards angled




JP2074RB	115 × 2.2,	反咬切, 左	left
JP2075RB	115 × 2.2,	反咬切, 右	right
JP2132RB	210 × 2.2,	反咬切, 左	left
JP2133RB	210 × 2.2,	反咬切, 右	right

注: 上述“左弯”、“右弯”为杆部左弯或右弯 12°; “上弯”为杆部上弯 15°; “上翘”为头部上翘 20°
Remark: Above left or right curved means shaft left or right curved to 12°; Upwards curved means shaft upwards curved to 15°; Upwards angled means tip upwards angled to 20°.

篮钳


Basket Forceps



- L × B**
- JP2076RB** 125×2.2, 侧咬切, 左 left
 - JP2077RB** 125×2.2, 侧咬切, 右 right
 - JP2078RB** 125×3.4, 侧咬切, 左 left
 - JP2079RB** 125×3.4, 侧咬切, 右 right
 - JP2125RB** 215×2.2, 侧咬切, 左 left
 - JP2126RB** 215×2.2, 侧咬切, 右 right

抓钳

Grasping Forceps




- JP2157RB** 125, 带锁 with lock
- JP2158RB** 215, 带锁 with lock
- JP2159RB** 85, 带锁 with lock



- JP2236RB** 125, 常规





- JP2235RB** 125, 常规
- JP2164RB** 215, 带锁 with lock



关节镜器械
 Arthroscopic instrument

关节钳
 Grasper





- JP2150RB** 125, 带锁 with lock
- JP2151RB** 215, 带锁 with lock
- JP2152RB** 90, 带锁 with lock



- JP2265RB** 155, 带锁 with lock

穿线钳
 Suture Forceps



- JP2242RB** 155, 直 straight
- JP2243RB** 155, 弯 (上弯 15°) upward curved 15°



- JP2244RB** 155, 弯 (左弯 45°) left curved 45°
- JP2245RB** 155, 弯 (右弯 45°) right curved 45°



- JP2246RB** 155, 弯 (上翘 20°) upward angled 20°
- JP2247RB** 155, 弯 (上翘 40°) upward angled 40°




- JP2248RB** 155, 弯 (左弯 45°, 下开口) left curved 45°, down opening
- JP2249RB** 155, 弯 (右弯 45°, 下开口) right curved 45°, down opening



关节镜器械
 Arthroscopic instrument



抓线器
Suture Grasper



JP2237RB 155

JP2259RB 155


剪线器
Suture Cutter



JP2252RB 155 (闭口)


JP2238RB 155 (开口)

推结器
Knot Pusher




JP2256RB 155

JP2257RB 220




关节镜器械
Arthroscopic instrument

关节剪
Scissors




JP2170RB 120, 指圈式

关节剪
Scissors



JP2175RB 120, 指圈式 (20°)
JP2178RB 120, 指圈式 (60°)

关节剪
Scissors




JP2173RB 125, 雪茄式 (20°)
JP2177RB 125, 雪茄式 (60°)



关节镜器械
Arthroscopic instrument


缝合钳
Suture Passers



骨牵引针
Guiding Wire




骨刀
Osteotome




- ZP2630RN** 160×4.5×15°, 前弯 forward curved
- ZP2631RN** 160×4.5×15°, 后弯 backward curved
- ZP2632RN** 160×4.5×30°, 前弯 forward curved
- ZP2633RN** 160×4.5×30°, 后弯 backward curved

骨锉
Bone File



- ZP2490RN** 165, 弯 (上) upward curved
- ZP2491RN** 165, 弯 (下) backward curved



关节镜器械
Arthroscopic instrument

结扎缝合引线器
Suture Guide



ZP2401RN 直型



ZP2404RN 左弯 (45°) left curved
ZP2405RN 右弯 (45°) right curved



ZP2406RN 左弯 (60°) left curved
ZP2407RN 右弯 (60°) right curved



关节镜器械
Arthroscopic instrument

骨科用穿线器
Suture Passer



ZP2410RN 155, 直 straight




ZP2415RN 155×45°, 左弯 left curved
ZP2416RN 155×45°, 右弯 right curved



ZP2419RN 145×90°, 左弯 left curved
ZP2420RN 145×90°, 右弯 right curved

骨科用穿线器
Suture Passer



- ZP2473RN 100×30°, 左弯 left curved
- ZP2476RN 100×30°, 右弯 right curved



- ZP2474RN 100×45°, 左弯 left curved
- ZP2477RN 100×45°, 右弯 right curved




- ZP2475RN 100×65°, 左弯 left curved
- ZP2478RN 100×65°, 右弯 right curved

穿孔针
Pin Punch




- ZP2471RN 80, 尖头 (直) straight
- ZP2472RN 80, 尖头 (弯) curved

小切口手术缝合线拉钩
Suture Catcher



- ZP2494RN 160, U型 U shaped
- ZP2493RN 220, 圆头 round tip



关节镜器械
Arthroscopic instrument



探针
 Probe



JP2180RB 115 (3)
JP2183RB 115 (5)

刮匙
 Curette



JP2213RB 115 (3)
JP2214RB 115 (5)

手术刀
 Knife



JP2191RB 110, 反角形

手术刀
 Knife



JP2193RB 110, 香蕉形

手术刀
 Knife



JP2194RB 110, 香蕉形

手术刀
 Knife



JP2196RB 110, 玫瑰形

手术刀
 Knife



JP2197RB 110, 镰刀形

手术刀
 Knife



JP2198RB 110, 钩形

手术刀
 Knife



JP2207RB 110, 半月板

手术刀
 Knife



JP2209RB 110, 韧带松弛

平锉
 File



JP2219RB 115 (细齿)

平锉
 File




JP2220RB 115 (粗齿)



关节镜器械
 Arthroscopic instrument

交叉韧带定位系统
Cruciate Ligament Positioning System



导向手柄
Handle



ZP2501RN 110

导向管
Guiding Sleeve




ZP2502RN $\Phi 3$

导向管
Guiding Sleeve



ZP2503RN $\Phi 2.5$

瞄准臂
Sighting Arm



ZP2510RN 点对点

瞄准臂
Sighting Arm



ZP2511RN 点对肘

瞄准臂
Sighting Arm




ZP2512RN 后叉胫骨

瞄准臂
Sighting Arm



ZP2513RN 后叉股骨



关节镜器械
Arthroscopic instrument

交叉韧带定位系统
 Craciate Ligament Positioning System



导向手柄
 Handle



ZP2700RN 110

导向管
 Guiding Sleeve




ZP2711RN Φ 2.5

导向管
 Guiding Sleeve




ZP2713RN Φ 3.5

瞄准臂
 Sighting Arm



ZP2721RN PCL 胫骨

瞄准臂
 Sighting Arm



ZP2722RN PCL 股骨

瞄准臂
 Sighting Arm



ZP2723RN ACL 股骨

瞄准臂
 Sighting Arm




ZP2724RN ACL 胫骨



关节镜器械
 Arthroscopic instrument

直径尺
Diameter Ruler




ZP2540RN 200, 测量器

后又保护器
Back Ligament Proteetor




ZP2563RN 120

测探器
Depth Gauge




ZP2561RN 240

取腱器
Tendon Extractor




ZP2572RN Φ 7



关节镜器械
Arthroscopic instrument

空心钻


Twist Drill for tibia



- ZP2590RN $\Phi 5$, 胫骨
- ZP2591RN $\Phi 6$, 胫骨
- ZP2592RN $\Phi 7$, 胫骨
- ZP2593RN $\Phi 8$, 胫骨
- ZP2594RN $\Phi 9$, 胫骨
- ZP2595RN $\Phi 10$, 胫骨

空心钻


Twist Drill for femur



- ZP2610RN $\Phi 5$, 股骨
- ZP2611RN $\Phi 6$, 股骨
- ZP2612RN $\Phi 7$, 股骨
- ZP2613RN $\Phi 8$, 股骨
- ZP2614RN $\Phi 9$, 股骨
- ZP2615RN $\Phi 10$, 股骨
- ZP2618RN $\Phi 4.5$, 股骨

骨钻

Bone Drill



- ZP2703RB $\Phi 7$
- ZP2704RB $\Phi 8$
- ZP2705RB $\Phi 9$




关节镜器械
 Arthroscopic instrument

股骨瞄准器
Femur Aimer




- ZP2522RN 5
- ZP2523RN 6
- ZP2524RN 7

定位针
Positioning Pin



- ZP2535RN $\Phi 2.4$
- ZP2530RN $\Phi 1$

骨牵引针
Guiding Pin




- ZP2280RN $\phi 4.3 \times 330$, 圆型
- ZP2281RN $\phi 4 \times 330$, 圆型 (中空) cannulated
- ZP2282RN $\phi 4.2 \times 330$, 圆型

显微刀
Micro Knife



- ZP2290RN 240×4.4 , 弯型 (可拆) detachable




关节镜器械
Arthroscopic instrument

微骨折锥
 Microfracture Awl



	JP2225RB	125, 直	straight	JP2230RB	215, 直	straight
	JP2226RB	125, 弯 (25°)	curved	JP2231RB	215, 弯 (25°)	curved
	JP2227RB	125, 弯 (40°)	curved	JP2232RB	215, 弯 (40°)	curved
	JP2228RB	125, 弯 (60°)	curved	JP2233RB	215, 弯 (60°)	curved
	JP2229RB	125, 弯 (90°)	curved	JP2234RB	215, 弯 (90°)	curved

微骨折锥
 Microfracture Awl



	ZP2661RN	120, 弯 (25°)	curved
	ZP2662RN	120, 弯 (40°)	curved
	ZP2663RN	120, 弯 (60°)	curved
	ZP2664RN	120, 弯 (90°)	curved

打入器
 Impactor



	ZP2670RN	100
--	-----------------	-----




骨刀
Osteotome



ZP2646RN 110×4, 弯 curved


ZP2647RN 110×6, 弯 curved

骨凿
Chisel




- ZP2651RN 240×4, 直 straight
- ZP2652RN 240×4, 弯 curved
- ZP2653RN 240×6, 直 straight
- ZP2654RN 240×6, 弯 curved

骨刮匙
Bone Curette Cannulated




- ZP2671RN 120×3.5, 空心
- ZP2672RN 120×4.5, 空心
- ZP2673RN 120×5.5, 空心

骨刮匙
Bone Curette Cannulated




- ZP2641RN 110×3.5×15°, 空心
- ZP2642RN 110×5×15°, 空心




关节镜器械
Arthroscopic instrument

打孔器
Punch




- ZP2742RN 120×Φ3.5
- ZP2744RN 120×Φ5

骨用丝锥
Screw Tap



- ZP2746RN Φ4.6×140
- ZP2747RN Φ4.8×140
- ZP2748RN Φ5.2×140

开口锥
Screwdriver



- ZP2751RN 140×2.7
- ZP2752RN 140×3.8
- ZP2753RN 140×5.8
- ZP2754RN 230×5



关节镜器械
Arthroscopic instrument

足踝微骨折器械

Ankle Microfracture Instruments



骨锥
Bone Awl



ZP2760RN 170, 弯 (25°) curved

骨锥
Bone Awl



ZP2761RN 170, 弯 (40°) curved

骨锥
Bone Awl



ZP2762RN 170, 弯 (60°) curved

骨锥
Bone Awl



ZP2763RN 170, 弯 (90°) curved

打入器
Impactor



ZP2764RN 170

骨刮匙
Bone Curette



ZP2765RN 170×5×15°, 双刃
ZP2766RN 170×3×15°, 双刃
ZP2767RN 170×2×15°, 双刃

骨刮匙
Bone Curette



ZP2768RN 170×5, 单刃
ZP2769RN 170×3, 单刃
ZP2770RN 170×2, 单刃

骨刀
Osteotome




ZP2771RN 170×15° (直) straight
ZP2772RN 170×15° (上弯) left curved
ZP2773RN 170×15° (下弯) downward curved
ZP2774RN 170×15° (窄形) narrow

探针
Probe



ZP2775RN 170×5



关节镜器械
Arthroscopic instrument



ZP2000 膝关节手术器械
 Knee joint instruments set

序号	产品编码	产品名称	规格型号	数量	单位
1	JP2003RB	篮钳	125×3.2, 鸭嘴, 直	1	把
2	JP2021RB	篮钳	125×3.2, 鸭嘴, 左弯	1	把
3	JP2022RB	篮钳	125×3.2, 鸭嘴, 右弯	1	把
4	JP2023RB	篮钳	125×3.2, 鸭嘴, 上弯	1	把
5	JP2007RB	篮钳	125×3.2, 鸭嘴, 上翘	1	把
6	JP2038RB	篮钳	125×4.4, 卵圆形	1	把
7	JP2075RB	篮钳	115×2.2, 反咬切	1	把
8	JP2074RB	篮钳	115×2.2, 反咬切	1	把
9	JP2076RB	篮钳	125×2.2, 侧咬切	1	把
10	JP2077RB	篮钳	125×2.2, 侧咬切	1	把
11	JP2150RB	关节钳	125, 带锁	1	把
12	JP2180RB	探针	115	1	把
13	JP2173RB	关节剪	125, 雪茄式	1	把
14	JP2198RB	手术刀	110, 钩形	1	把



关节镜器械
 Arthroscopic instrument



ZP2100 肩关节手术器械
 Shoulder joint instruments set

序号	产品编码	产品名称	规格型号	数量	单位
1	JP2242RB	穿线钳	155, 直	1	把
2	JP2244RB	穿线钳	155, 弯	1	把
3	JP2245RB	穿线钳	155, 弯	1	把
4	JP2247RB	穿线钳	155, 弯	1	把
5	JP2237RB	抓线器	155	1	把
6	JP2252RB	剪线器	155	1	把
7	JP2256RB	推结器	155	1	把
8	JP2265RB	关节钳	155, 带锁	1	把
9	ZP2490RN	骨锉	165, 弯	1	把
10	ZP2631RN	骨刀	160×4.5×15°, 后弯	1	把
11	ZP2494RN	小切口手术缝合线拉钩	160, U型	1	把
12	JP2183RB	探针	115	1	把
13	ZP2280RN	骨牵引针	Φ4.3×330, 圆型	1	把
14	ZP2401RN	结扎缝合引线器	直型	1	把
15	ZP2404RN	结扎缝合引线器	左弯	1	把
16	ZP2405RN	结扎缝合引线器	右弯	1	把
17	ZP2406RN	结扎缝合引线器	左弯	1	把
18	ZP2407RN	结扎缝合引线器	右弯	1	把



关节镜器械
 Arthroscopic instrument




ZP2200 交叉韧带重建器械
Cruciate ligament Reconstruction instruments set

序号	产品编码	产品名称	规格型号	数量	单位
1	ZP2501RN	导向手柄 DXSB-110	110	1	把
2	ZP2502RN	导向管 DXG-3	Φ3	1	把
3	ZP2503RN	导向管 DXG-2.5	Φ2.5	1	把
4	ZP2510RN	瞄准臂 MZB-D	点对点	1	把
5	ZP2511RN	瞄准臂 MZB-Z	点对肘	1	把
6	ZP2512RN	瞄准臂 MZB-JG	后叉胫骨	1	把
7	ZP2513RN	瞄准臂 MZB-GG	后叉股骨	1	把
8	ZP2522RN	股骨瞄准器 GMQ-5	5	1	把
9	ZP2523RN	股骨瞄准器 GMQ-6	6	1	把
10	ZP2540RN	直径尺	200, 测量器	1	个
11	ZP2561RN	测深器 CSQ-240	240	1	把
12	ZP2563RN	后叉保护器 BHQ-120	120	1	把
13	ZP2572RN	取腱器 QJQ-7	Φ7	1	把
14	ZP2493RN	小切口手术缝合线拉钩	220, 圆头	1	把
15	ZP2618RN	空心钻 GGZ-4.5	Φ4.5, 股骨	1	个
16	ZP2611RN	空心钻 GGZ-6	Φ6, 股骨	1	个
17	ZP2612RN	空心钻 GGZ-7	Φ7, 股骨	1	个
18	ZP2613RN	空心钻 GGZ-8	Φ8, 股骨	1	个
19	ZP2614RN	空心钻 GGZ-9	Φ9, 股骨	1	个
20	ZP2615RN	空心钻 GGZ-10	Φ10, 股骨	1	个
21	ZP2591RN	空心钻 JGZ-6	Φ6, 胫骨	1	个
22	ZP2592RN	空心钻 JGZ-7	Φ7, 胫骨	1	个
23	ZP2593RN	空心钻 JGZ-8	Φ8, 胫骨	1	个
24	ZP2594RN	空心钻 JGZ-9	Φ9, 胫骨	1	个
25	ZP2595RN	空心钻 JGZ-10	Φ10, 胫骨	1	个
26	ZP2535RN	定位针 DWZ-2.4	Φ2.4	1	个
27	ZP2530RN	骨牵引针	Φ1	1	个




关节镜器械
Arthroscopic instrument



ZP2400 微骨折手术器械
Microfracture instruments set

序号	产品编码	产品名称	规格型号	数量	单位
1	ZP2760RN	骨锥	170, 弯	1	把
2	ZP2761RN	骨锥	170, 弯	1	把
3	ZP2762RN	骨锥	170, 弯	1	把
4	ZP2763RN	骨锥	170, 弯	1	把
5	ZP2764RN	打入器	170	1	把
6	ZP2765RN	骨刮匙	170×5×15°, 双刃	1	把
7	ZP2766RN	骨刮匙	170×3×15°, 双刃	1	把
8	ZP2767RN	骨刮匙	170×2×15°, 双刃	1	把
9	ZP2768RN	骨刮匙	170×5, 单刃	1	把
10	ZP2769RN	骨刮匙	170×3, 单刃	1	把
11	ZP2770RN	骨刮匙	170×2, 单刃	1	把
12	ZP2771RN	骨刀	170×15°	1	把
13	ZP2772RN	骨刀	170×15°	1	把
14	ZP2773RN	骨刀	170×15°	1	把
15	ZP2774RN	骨刀	170×15°	1	把
16	ZP2775RN	探针	170×5	1	把



关节镜器械
Arthroscopic instrument

