

SAS  
series

## Ampoule Steam Sterilizer


SAS系列 安瓿蒸汽灭菌器



### Product Introduction 产品简介

SAS安瓿蒸汽灭菌器采用高温高压的纯蒸汽作为灭菌介质，可实现101°C~130°C温度区间的均匀灭菌，广泛用于安瓿、口服液及小容量静脉输液等药品制剂的灭菌、检漏和清洗处理。

SAS ampoule steam sterilizer adopts hot and pressurized pure steam as the sterilizing medium to achieve the uniform sterilization between 101~130°C. It is widely used for sterilization, leak detection and washing of ampoules, oral liquid and small volume IV products.



安瓿蒸汽灭菌器

### Product Advantages 产品优势

- 标准化的管路系统，结构紧凑，便于维护；  
Standardized pipeline system, compact structure, easy maintenance;
- 多重检漏方式，检漏更彻底；全新的清洗方式，清洗更高效节能；  
Multiple leak detection methods make leak detection more thorough; brand new washing methods make washing more efficient and energy-saving;
- 机器人自动焊接、磨光工艺，质量更可靠；  
Automatic robot welding, automatic polishing, reliable quality;
- 格栅分层式装载内车，温度均匀性更高；  
Sterilization trolley with layered rack makes temperature distribution more uniform;
- 专用控制系统，控温更精确。  
Dedicated control system for more accurate temperature control.



### Product Use 产品用途

- 玻璃安瓿注射剂  
Glass Ampoule Injection
- 西林瓶注射剂  
Vial Injection
- 玻璃瓶注射剂  
Glass Bottle Injection

